



Penn Medicine



Center for
Digital Health

Digital Data and Mental Health Study

Electronic Communication and Digital Data

Downloading and Donating Data

Hello! This document provides an overview of how to safely and securely download and donate your data

You can donate independently or work with a CDH researcher.

If you have any questions or concerns, please email

digitalhealth@pennmedicine.upenn.edu



You will download AWARE Client V2 app on your Android

AWARE Client V2 securely and passively collects screen time activity and movement from your Android



How to Download AWARE Client V2

Google's recent move to block accessibility services from being published in the Play Store has led to a temporary removal of AWARE's main client and plugins from the Play Store.

The following is a workaround to install AWARE if you search for it in the Google Play store and it cannot be found.

Step 1 Configure Phone for Installation of Apps Outside the Play Store

Go to "Settings -> Apps & notifications -> Advanced -> Special app access -> Install unknown apps"

Select the app (your preferred internet browser) that you wish to grant permissions to install from unknown sources from.



How to Download AWARE Client V2

Step 2: Click on the link below from your phone to download the AWARE V2 app


<http://jenkins.awareframework.com/job/com.aware.phone/lastSuccessfulBuild/artifact/aware-phone/build/outputs/apk/release/aware-phone-armeabi-release.apk>



How to Download AWARE Client V2

Step 3: Open and install the AWARE V2 App


1. Install AWARE V2 by clicking “INSTALL”
2. Open the AWARE V2 app after installation is complete
3. Enable all permissions (press “Allow”) as the messages appear






How to Download AWARE Client V2



Step 4 You will see this screen. Go to Settings



Step 5 Click the QR code symbol



Step 6 Click 'OK'





How to Download AWARE Client V2

Step 7 Take a picture of this QR code

The QR code is a unique code for this study and pre-populates AWARE's data sources.





How to Download AWARE Client V2

Step 8 Select 'Join'





Step 9: There may be more permissions requested depending on the study. Press "Allow" for each.



How to Download AWARE Client V2



Step 10: Under 'Settings You will see this screen, please click 'Device Name'





How to Download AWARE Client V2



Step 11 Please change 'Device Name' to 'Your Unique ID Number' from the Center for Digital Health Page





How to Download AWARE Client V2


Step 12 Click refresh to update the Device's Name. Please click 'Close'






How to Download AWARE Client V2

Step 13 Click 'Advanced Settings'



Step 14, Select Off for WiFi, Battery Charging Only, and select "Hide settings" for UI Mode




*Note, once you return back to **"Settings"** you will no longer be able to view or access **"Advanced Settings"** from your phone.





How to Download AWARE Client V2

Step 15 Click “Settings” to return back to its original page, and then click the Cloud icon in the top left part of the screen to sync your data.




Step 16 Read the following information, and then click “Execute.” You will see this screen when completed, please click ‘close’



Congrats!! You're finished!!! You've donated:



AWARE Client V2



[This Photo](#) by Unknown Author is licensed under
[CC BY-NC-ND](#)



Penn Medicine



Center for
Digital Health

**THANK
YOU!**

If you have any questions or concerns please email



digitalhealth@penntmedicine.upenn.edu